


Getting Started on Advisor Pro and DreaMed Diary app

Setting Up Your Advisor Pro Account

- 1 Your healthcare professional will send you a welcome email to Advisor Pro. Access the email from your phone or pc and tap on the link to set up your account password.
Note:
 - The email link is valid for 72 hours. If it expires, your healthcare professional must send you a new one.
 - Password must be at least 8 characters, contain at least one upper-case letter, lower-case letter, and a number or a special character.

- 2 Confirm your acceptance of the Advisor Pro terms of use and Privacy Policy.

- 3 Enter your credentials and tap on **Sign Up**



Confirmed Email:
[xxxx@gmail.com](#)

Please select your new password:

Password

Confirm Password

I've read and agreed to Advisor Pro
[Privacy Policy](#) and [Terms of Use](#)

Sign Up

- 4 Download DreaMed Diary app to your smartphone



- 5 Open the app by tapping on the icon and login using your email and password




Dairy App Use and Connecting Dexcom Account

1


The following screen will be displayed upon first launching the app. Tap on the green button to **Connect Your Dexcom Account**

Note: This can also be done at a later time through the settings menu under **Connected Devices**.



2

Enter your Dexcom account username and password and **Login**



3

Type in your full name and then tap the **Authorize** button to allow DreaMed Diary app to automatically upload information from your Dexcom CGM into your diary app

At the bottom of your app screen you will find a navigation bar for 3 different menus, tap on the icons to navigate between screens, red indicates the opened menu currently showing on your screen:




Plan Settings

The Plan Settings is configured by your healthcare provider. You can view your plan here.

Tap on the different periods of the day to view the bolus plan for the **Morning, Afternoon, Evening, and Night.**

From time to time, your HCP may decide to modify your treatment plan using the Advisor Pro. You will receive notification that there is a new plan awaiting your review and approval. Only after you accept the new plan, it will be updated in your plan settings.

For further guidance on reviewing and changing your treatment plan, please refer to the guide “Understand Your Treatment Plan” in the educational portal.



Plan Settings Edit

Basal plan (long acting insulin)

10:30 PM 19 Units
Basaglar

Bolus plan (short / rapid acting insulin)

Morning **Afternoon** Evening Night

Afternoon **11 AM — 5 PM**

Insulin type Novolog

Carbs ratio 10 g/u


Correction factor 15 mg/dL/u

Glucose target 115 mg/dL

Diary Plan Settings

Settings Menu

The DreaMed Diary app includes various settings that define its behavior and how some things appear on the screen, such as how carb events are reported. Typically, your HCP will set these things up for you. If needed, you can modify them in the **Settings** screen. However if you are unsure, discuss with your healthcare provider before making changes to your plan.



Settings

Profile →

General Settings →

Connected Devices →

About →

Help →

Diary Plan **Settings**


In the General Settings, if you count carbs per meal, you can modify your bolus calculator settings which is made available by your health care provider.

Connect your Dexcom CGM using “add new device” and following the instructions at the beginning of this guide.

Diary Screen

The Diary log screen displays a chronological list of all the events. Each event appears with its specific date and time.

Note: The “calculate dose” button will be available only if your health care provider has opened this option.



Adding events to your diary

Note: Events may be entered manually as described below.

- 1 To add an event, go to the diary screen and click on **Add Event**
- 2 In the opened menu enter the information you would like to capture, as described in the picture below. You do not have to enter all of the fields in every event. Enter only the relevant information to the specified date and time.

For recording insulin injections tap the + button, select the type of insulin and enter the number of units injected.

Record your meal size or carbs amount consumed in your meal depending on your treatment type.

Enter any additional information you would like to record for this event.

The screenshot shows the 'Add Event' form with the following fields and callouts:

- Date and Time:** 'Today, May 13' and '08:52'. A callout explains that this information will be assigned to the new diary event and can be changed by tapping the row.
- Insulin:** A field with a purple insulin icon, the text 'Insulin', and a '+ Add another' button. A callout explains that users should tap the '+' button to select insulin type and enter units.
- Glucose Level:** A field with a red glucose icon, the text 'Glucose Level', a numeric input '0', and the unit 'mg/dL'. A callout explains that users can capture glucose readings as they appear on their glucometer device.
- Meal size:** A field with a blue fork and knife icon, the text 'Meal size', and a dropdown menu set to 'None'. A callout explains that users should record meal size or carb amount.
- Activity:** A field with a green running icon, the text 'Activity', and a dropdown menu set to 'None'. A callout explains that users should enter the intensity of the activity.
- Note:** A field with a document icon and the text 'Tap here to add a note'. A callout explains that users should enter any additional information here.

The date and time shown here will be assigned to the new diary event you are adding. To change it, simply tap on the row and chose the correct date and time from the display.

Here you can capture your glucose readings as they appear on your glucometer device.

Enter the intensity of the activity you performed relative to your usual physical activity.

- 3 When completed, click on **Save** at the top right corner to approve. This event now appears in your **Diary** screen

Edit/ Delete events

To edit or delete an event from your diary log, tap on the relevant event in your Diary screen, then tap on the **Edit** at the top right corner.

To edit: Change the event details as needed and tap on the **Save** at the right corner to approve.

To erase: Simply tap **Delete**

Note: Deleted events cannot be recovered

To return without editing/deleting tap on the **X** at the top left corner.

X **Event** **Delete** **Save**

Today, Feb 11 14:08

+ Insulin **+** Add another

x Novolog 20 Units

+ Glucose Level 0 mg/dL

+ Carbs Amount 0 Grams

+ Activity None

☰ Tap here to add a note